

Survey Of Bulk Tank Milk From All Danish Dairy Herds In 2009 And 2010 With Real time PCR







Torben W. Bennedsgaard



Foulum, Aarhus University





Our Milk

- a pure pleasure



Denmark - The world's number 1





Denmark - The world's number 1



Caroline Wozniacki - Miss Sunshine



Highly industrialised Dairy Industry







Number of dairy cows 573,000

Herd size 132 cows Yield 9,079 kg









The World's leader in robotic milking

AMS herds 22 %

Cows 27%










Organic herds 10%

10 % of milk





Herds in the Danish B-register from 1954 to 2010 in percent of all herds and from 2005 of active dairy herds





Aim of study

- To provide a sensitivity comparison of bacterial culture and the PCR assay (Pathoproof®) in Streptococcus agalactiae detection in Bulk tank milk samples from all 4258 Danish Dairy herds in 2009.
- To validate the PCR testing of all 12 genes as part of the annual nationwide BTM screening program

Bulk tank milk



Culture
Streptococcus agalactiae + / -





Collection and sampling



PCR – real time Ct value for 12 gene








Table 1. Culture and PCR for *Streptococcus agalactiae* from all 4258 Danish dairy herds in 2009

Test for GBS	PCR		
Culture	Positive	Negative	
Positive	178	20	
Negative	132	3928	

Table 2: Number of dairy herds PCR and culture positive for *Streptococcus agalactiae* 2009 and 2010

Year	herds	PCR GBS-positive	Culture GBS-Positive
2009	4258	301 (7.3 %)	198 (4.7 %)
2010	4093	271 (6.6 %)	141 (3.4 %)

Table 3: Results of PCR test PathoProof® in BTM samples from 4258 herds in 2009

Bacteria/gen	% NoCt	Lowest	10% Percentile
Staph aureus	9	20	29
Staph. Sp	0	18	27
Beta-lactam	22	22	31
Str. agalactiae	93	17	26
Str. dysgalactiea	14	16	28
Str. uberis	5	14	26
Coryne. bovis	10	25	32
Enterococcus	22	21	30
E.coli	39	18	30
Klebsiella	87	19	31
S. marcescens	98	25	34
A.pyo./ P. ind	37	19	32

Conclusion

- Streptococcus agalactia is an increasing problem in Denmark
- PCR is more sensitive than culture
 - Sensitivity 94%
- PCR on BTM samples can be used as surveiliance of other mastitis pathogens in dairy herds, and give farmers good information on further prophylactic actions against udder infections